MICROTEST SR

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
Fort

MICROTEST SR

Detection of anaerobic sulphate-reducing bacteria

DESCRIPTION

Each packing unit consists of 18 tubes with the SRB-medium and 2 tubes containing the respective number of sterilised capillaries. The SRB-medium may be stored at room temperature if protected from the light. It has a shelf life of 10 months.

PROCEDURE

- 1- Open the culture tube and stand it upright;
- 2- Using forceps, take one capillary from the container (if possible, flame the jaws of the forceps before use);
- 3- Dip the capillary into the fluid to be tested;
- 4- Using an absorbent paper, empty the capillary;
- 5- Insert capillary into the culture medium. Reseal the tube tightly;
N.B.: handle the capillary only with forceps. Do not touch it with your hands or clothing and do not lay it down. Use a fresh capillary for each sample;
- 6- Mark the necessary information on the tube and incubate at 35-37°C;
- 7- Compare the degree of blackening with the model chart below, at regular intervals, for 5 days;
- 8- A strong contamination will be evident after approximately 48 hours; a slight contamination will reveal after 72 hours;
- 9- For comparative testing, development time and temperature should be the same.

INTERPRETATION

The growth of anaerobic sulphide generating bacteria results in formation of black iron sulphide. Blackening may begin at any point of the tube, its intensity depending on the degree of contamination. The result is not quantitative but shows the approximate level of contamination.